

# Cozy™ 2.0

**Far Infrared Throw/Wrap**



**The New Far Infrared Cozy™ 2.0 wrap/throw uses  
Our latest most advanced Micron 2.0 Far infrared Technology.**

## Easy to Use

The Cozy 2.0 is designed to use while you sit relaxing, reading, watching TV or when you are just having a nap.

Our Cozy 2.0 is so versatile you can cover yourself and take your Far infrared treatment while you sit or lie comfortably on the bed, chair, sofa, treatment couch, while relaxing or working.

## Relax

Stress is a major issue for most of us and relaxation and de-stressing is essential for our lasting wellbeing and good health.

Your CozY™ gently warms the entire body creating a feeling of deep relaxation.

- Uses our latest Micron 2.0 Far Infrared heating system
- Super Lightweight
- Totally Flexible
- Ultra Low EMF's.
- †Low Odour – Good for Chemical Sensitivity
- Easy to clean Polyester
- †Temperature Range: 40-55°C
- Warranty – 12 Months

## Finest Quality Components

The Far Infrared CozY™ is unique and uses our latest, Micron 2.0 Far Infrared Heating technology.

## What You Get

**Controller & Leads:** adjusts and controls the temperature of the Far Infrared CozY™ 2.0. This connects directly to the mains for power via an electrical lead.

Telephone: 01727 602022

Website: [www.get-fitt.com](http://www.get-fitt.com)

Email: [admin@get-fitt.com](mailto:admin@get-fitt.com)

Select **Green** for 40-45°C, **Blue** for 45-50°C and **Red** for 50-55°C.



### Ending Your Session

When you wish to stop your session, just press and hold the **Power** button on the controller until it flashes and returns to ready mode; the CozyY™ 2.0 will now stop heating, then **unplug the mains lead**.

### Storage

Your Far Infrared CozyY™ must be stored in a dry and cool place.

### What You Get

Controller: Adjusts and controls the temperature and timer settings.

Universal Adaptor: Connects the Controller to the Mains.

Mains Leads: Connect directly to the mains for power via an electrical lead.

## Components

### AC/DC Universal Adapter



## Features

Micron 2.0 Far Infrared system

Super Lightweight

Totally Flexible

Ultra-Low EMF's.

Low Odour – Good for Chemical Sensitivity

Easy to clean Polyester

Temperature Range: 40-55°C.

Warranty – 12 Months

The Far Infrared CozY™ 2.0 is unique and uses our most advanced Ultra Low EMF Micron 2.0 Far Infrared technology.

Telephone: 01727 602022

Website: [www.get-fitt.com](http://www.get-fitt.com)

Email: [admin@get-fitt.com](mailto:admin@get-fitt.com)

## Specifications

- Colour: Grey
- Material: Polyester
- Input Volts: 100V ~240V 50/60Hz 1.5A
- Output: DC 24V 2.0A
- Output Power: 48W
- **Control Panel:** Manual Control
- Temperature: 40°C - 55°C
- Dimensions: 91 x 150 X 0.5 cm
- Net Weight: 1.0 Kg
- Auto Shut Off: 4 Hours
- **Far Infrared Emissions:** 5-20 microns
- **Far Infrared Emitter:** Get Fitt – Micron 2.0 Technology
- Warranty 12 Months

Telephone: 01727 602022

Website: [www.get-fitt.com](http://www.get-fitt.com)

Email: [admin@get-fitt.com](mailto:admin@get-fitt.com)