

# COZY 4.0

## Far Infrared Throw/Wrap



**The Far Infrared Cozy™ 4.0 Wrap and Throw uses  
Our *latest* most advanced Micron 4.0 Far infrared Technology.**

# Far Infrared COZY 4.0

## Safety Essentials

**Do not Fold while in use.**

**Do not use in the BED & sleep on the Cozy 4.0.**

**Unplug the power cord from the mains after each use.**

### **Easy to Use**

The versatile Cozy 4.0 is designed to be used to cover you while you relax sitting reading, watching TV or when you are just having a nap.

### **Relax and Detox**

The deep penetration of the Far Infrared rays offer an extremely comfortable and efficient means of body cleansing and detoxification

Your Cozy 4.0 gently warms the entire body creating a feeling of deep relaxation.

### **Storage**

Your Far Infrared Cozy 4.0 must be stored in a dry and cool place.

### **Finest Quality Components**

The Far Infrared Cozy 4.0 uses our *latest* most advanced Micron 4.0 Far infrared Technology.

### **What You Get**

Cozy 4.0 with Controller: Adjusts and controls the temperature.

## Using Your Cozy 4.0

### 1. Plug the connector into the Cozy 4.0.socket: 3 Steps

Step 1



Step 2



Step 3



**2. Plug the controller mains cable with the 3 pin plug into the mains and then switch on the power at the mains plug. The display light will come on - The Cozy 4.0 is **not** heating yet.**

Temperature Display



Power Button & Temperature control



Telephone: 01727 602022

Website: [www.get-fitt.com](http://www.get-fitt.com)

Email: [admin@get-fitt.com](mailto:admin@get-fitt.com)

### 3. Turning on Your Cozy 4.0

Press and hold the selector button until a number shows on the display; then release the button.

**Press & Hold the Power Button**



### 5. Press the selector to choose your temperature.

**Press the Temperature Button to select your temperature**



### Temperature Settings

Setting	1	2	3	4	5	6	H
Temperature °C	35	40	45	50	55	60	65

Once activated your Cozy 4.0 can take 2/3 minutes to reach the desired temperatures, this can take a little longer for the higher operating temperatures and may vary according to the ambient room temperature.

## **Ending Your Session**

When you wish to stop using the Cozy 4.0, press and hold the power button until the display shows a horizontal line. Then unplug when finished.

**NOW unplug from the mains**

## Washing Instructions

The Cozy 4.0 can be gently hand washed in warm water, or machine washed at a maximum temperature of no more than 20 degrees Celsius; ideally using a non-biological liquid soap

It should only be HANG DRIED AND NOT MACHINE DRIED.

### Features

Micron 4.0 Far Infrared system

Super Lightweight

Totally Flexible

Ultra-Low EMF's.

Easy to clean Polyester

Temperature Range: 30-65°C.

Warranty – 12 Months

The Far Infrared Cozy™ 4.0 is unique and uses our most advanced Ultra Low EMF Micron 4.0 Far Infrared technology.

Telephone: 01727 602022

Website: [www.get-fitt.com](http://www.get-fitt.com)

Email: [admin@get-fitt.com](mailto:admin@get-fitt.com)

**Safety Instructions– Warning (Please read carefully)**

Do NOT use while folded.

Unplug the power cord after each use.

When sitting or lying on the unit use on Green/ Blue Temperature Settings.

DO NOT use in the bed and sleep on the CozY 4.0.

Do not use in or around water for danger of electrical shock.

Keep the electrical cord away from heaters or radiators.

Never operate the unit if the cord or plug has been damaged.

Do not pull on the electrical cord.

An authorized Dealer must perform any repair.

Check with your doctor regarding the use of this product.

Never cover the controller.

This product is not for medical purposes and is only manufactured for general use.

## Specifications

- Colour: Grey
- Material: Polyester
- Input Volts: 220V ~240V 50/60Hz 1.5A
- Output: DC 24V 4.0A
- Output Power: 400W
- **Control Panel:** Manual Control
- Temperature: 30°C - 65°C
- Dimensions: 180 x 80
- Net Weight: 1.0 Kg
- **Far Infrared Emissions:** 4-20 microns
- **Far Infrared Emitter:** Get Fitt – Micron 4.0 Technology
- Warranty 12 Months

Telephone: 01727 602022

Website: [www.get-fitt.com](http://www.get-fitt.com)

Email: [admin@get-fitt.com](mailto:admin@get-fitt.com)